

ERMA


MICROTOMES

ERMA MICROTOMES/DISPOSABLE MICROTOME BLADES



ERMA INC.
TOKYO JAPAN

SLEDGE MICROTOMES



Large Sledge Microtome : ESM-100L

Slide way length : 40cm

Maximum opening of object clamp : 39 × 43mm

Range of graduations for cutting of specimen : 0.5~12 μm in increments of 0.5 μm



Small Sledge Microtome : ESM-150S

ESM-100L / ESM-150S

Combination (selectable)
 Microtome Knife 17cm1
 Wooden handle for knife1
 Dispo.Blade Holder1

Slide way length : 30cm
 Maximum opening of object
 clamp : 27 × 32mm
 Range of graduations for
 cutting of specimen : 2~30nm

Microtomes

MODEL	DESCRIPTION	CAT.NO.
ERM-200P	Precision Rotary Microtome	08-200-0
ERM-230L	Large Rotary Microtome	08-230-0
ESM-100L	Large Sledge Microtome	08-100-0
ESM-150S	Small Sledge Microtome	08-150-0
EFM-300*	Freezing Microtome	08-300-0
EFM-500	Cylinder Microtome	08-500-0

Disposable Microtome Blades



ERMA's PATHO CUTTER is superior disposable blades with essential sharpened and special resin coated edge for exceedingly slicing various hardness of paraffined tissue blocks.

Disposable Knife Holder



140mm holder is mainly for Rotary Microtome use, and 220mm holder is for Sledge Microtome use.

Relevant Equipments

MODEL	DESCRIPTION	CAT. No.
F-1	Paraffin Section Heater	08-511-0
FB-2	Tissue Floating Bath, 2 baths	08-513-0
FB-1	Tissue Floating Bath	08-514-0
	Carbon Dioxide Freezing Apparatus for Sledge type	08-520-0
	Carbon Dioxide Freezing Apparatus for Freezing type	08-525-0
	Disposable Knife Holder for Rotary type	08-540-0
	Disposable Knife Holder for Sledge type	08-545-0
Patho Cutter 22°	Disposable Microtome Blades 22°	08-622-0
Patho Cutter I 35°	Disposable Microtome Blades 35°	08-635-0
Patho Cutter II 35°	Disposable Microtome Blades 35°	08-635-1

* Specifications subject to change without notice.



ERMA INC.
TOKYO JAPAN