

# CRP Ultra TurbiPak


Giving benefits of hs-CRP at the price of regular CRP

CRP is one of the most sensitive indicators of inflammation. Measurement of CRP is useful for the detection and treatment monitoring of infection, tissue injury, inflammatory disorders and associated diseases.

### SIGNIFICANCE OF hs-CRP (HIGH SENSITIVE CRP) MEASUREMENT

Detection of CRP levels below 3 mg/L in normal individuals gives us valuable information in terms of the risk of the future cardiovascular events.

- Low Risk < 1 mg/L
- Average Risk 1.0 to 3.0 mg/L
- High Risk > 3.0 mg/L



The CRP Ultra TurbiPak Reagent Kit is Used for the quantitative measurement of C-reactive protein by immunoturbidimetry method.

## Pack size



### CRP Ultra TurbiPak:

Product code : AGD-CRP50

Packsize: R1: 2 X 12.5 mL

R2: 2 X 12.5mL



### CRP Ultra Calibrator:

Product code : AGD-CRPCAL2.5

Packsize: 5 X 0.5mL

## Features of CRP Ultra TurbiPak :

- Ready to use liquid stable reagent
- Wide linear range: 0.1 mg/L to 420 mg/L.
- Highly sensitive kit giving accurate and precise measurements at very low levels of CRP
- It helps in reliable assessment of potential cardiovascular risk of normal healthy individuals
- Calibrators standardized against ERM - DA 474/IFCC
- Excellent correlation with a nephelometric assay
- Economical: Providing significant cost savings
- No significant interference by lipids, bilirubin and hemoglobin
- No prozone phenomenon up to 1000 mg/L concentration
- Can be used in various automated analyzers

**AGD Biomedicals (P) Ltd.**  
Mehta Trade Center,  
Sir M.V. Road, Andheri East,  
Mumbai - 400 099 INDIA

P : +91-22-28231061/66  
+91-22-28257999  
sales@agdbio.com  
www@agdbio.com

ISO : 9001 : 2015  
ISO : 13485 : 2016

